

Aluminum Heat Sink
 X-axis = 4"
 Y-axis = 3"
 Z-axis = 3/4"
 e.g. McMaster Carr p/n
 8975K231 or similar

Important note: depending on the specific microscope stage and objective working distance, the block dimensions listed here will need to be reduced.

MICROSCOPE-MOUNTED THERMAL STIMULATOR
Dimensional drawing
Astra Bryant
7/13/2023
astrab@uw.edu

Channels on surface for peltier leads (24 AWC)
 Wire location is approximate - use device for final measurements

